

CE
2460

EYECRYLTM Sert

Advanced Monofocal IOL

Inspired by

natural vision

mono **MORE** *technology*

biotech

VISIONCARE

REDEFINING MONOFOCAL SEGMENT

- A new Advanced Monofocal IOL with Aspheric surface
- Continuous vision from distance to intermediate
- Spectacle independence for Intermediate Vision
- An extended range of vision: >1.75 D at spectacle plane with Visual Acuity 0.2 LogMAR or better
- Intermediate vision up to 57 cm
- MICS (2.2 mm) compatible to reduce SIA
- Designed for fast neural adaptation
- Minimal glare and haloes
- ABBE no. 47



PLHFD6



PLHFD6T


ADDED INTERMEDIATE VISION:

mono **MORE** *technology*


UNIQUE ASYMMETRIC POWER DISTRIBUTION

- Designed for extended vision
- Optimized optical Zones to provide extended depth of focus
- Asymmetric power distribution to minimize pupil dependency
- Photic phenomena similar to Monofocal IOLs

MONOFOCAL IOL




ADVANCED MONOFOCAL IOL



OPTIMIZED LIGHT ENERGY


Optimized optical zones provide the extended depth of focus for Intermediate vision and asymmetric power distribution minimizes pupil dependency in all lighting conditions.



DEFOCUS CURVE*


Extended Depth of focus for Intermediate vision without compromising far vision

- Continuous vision from Distance to Intermediate
- Great support for daily activities
- 0.2 lo MAR visual acuity even at +1.0D denotes tolerance range of IOL



MODULATION TRANSFER FUNCTION*


- Optimized energy distribution from Distance to Intermediate
- More than 0.3 MTF (50 lp/mm) @ 3.00 mm aperture
- Good contrast in all light conditions



OPTICAL RESULTS*

USAF IMAGES AT 3.0 MM APERTURE

Optimized light distribution to maintain better resolution and good contrast sensitivity



PRE-LOADED DELIVERY SYSTEM


SINGLE HAND IMPLANTATION WITH CONTROL OF SCREW TYPE INJECTOR:

BEST OF BOTH WORLDS



SIMPLE IOL IMPLANTATION USING FOLLOWING STEPS

STEP 01



Completely fill the entire length of the lens with BSS (don't use chilled BSS) until the entire lens is wrapped and lubricated. Also add lubrication to the tip of the injector until it is completely filled.

STEP 02



Push the blue injector plunger forward until the front push plate is flush against the injector housing. It is recommended to execute this step slowly and gently.

STEP 03



Inject adequate amount of any Biotech certified OVD having low to moderate viscosity, as shown here. Hydroxypropyl methylcellulose is recommended. The OVD should flow up to leading haptic of the IOL. Inject OVD from tip of the cartridge also, to fill the cartridge nozzle. Do not completely fill the chamber as this can move the IOL during insertion.

STEP 04




Close the cartridge flaps. Ensure that the flaps are locked with a "Click" sound.

STEP 05



Push the blue injector plunger forward until the rear push plate is flush against the injector housing or until the drive wheel of the injector moves. It is recommended to execute this step slowly and gently.

STEP 06



Hold the delivery system with a "Pen Grip", as shown here and keep your index finger on Drive Wheel.

STEP 07



Hold the system with the cartridge tip in a bevel down position. Now using your index finger, pull and rotate the drive wheel back slowly in order to push the lens forward until it is delivered.


MONOFOCAL IOL VS EYECRYLTM Sert IOL

Advanced Monofocal IOL

MTF COMPARISON*

Monofocal IOL provides 0.5 D to 0.8 D depth of field and EYECRYLTM Sert IOL provides more than 2.0 D depth of field (>1.75D depth of field at spectacle plane) which is helpful for Intermediate distance range of daily visual activities.


MONOFOCAL IOL



Modulation Transfer Function @3.0 mm aperture

EYECRYLTM Sert IOL

Advanced Monofocal IOL




Modulation Transfer Function @3.0 mm aperture

MONOFOCAL IOL VS. EYECRYL[™] Sert IOL

Advanced Monofocal IOL


USAF IMAGES COMPARISON*

MONOFOCAL IOL



EYECRYL[™] Sert IOL


Advanced Monofocal IOL



POST OPERATIVE OUTCOMES* (n=70)

SPHERICAL RESIDUAL


Post-operative spherical result shows that 97.14% of eyes were within ± 0.50 D spherical residual. All eyes were within ± 0.75 D spherical residual.




SPECTACLE INDEPENDENCE FOR DIFFERENT DISTANCES*

As per post-operative data, satisfaction for spectacle independence was found 100% for far vision & 91.89% for intermediate vision.

FAR VISION




INTERMEDIATE VISION



*Data on File.


PHOTIC PHENOMENON & VISUAL DISTURBANCE*

No case of dysphotopsia or post-operative disturbance related to vision have been reported so far. With minimal glare and halos, photic phenomenon is similar to monofocal IOLs.




PATIENT SATISFACTION*

Patient satisfaction was found to be very high. 98.88% patients would choose the same IOL again and recommend it to others.




REFRACTIVE CYLINDRICAL RESIDUAL

The cylindrical accuracy was within $\pm 0.5D$ for 85.71% cases. These post-operative cylindrical results prove the excellent corneal astigmatism correction.




IOL ROTATIONAL STABILITY

The lens has excellent rotational stability due to its advanced design and material. 92.59 % patients having the rotation within 5 degree proves the excellent rotational stability of the IOL.



QUALITY OF VISION*

100% patients were satisfied for their far prominent day to day activities e.g. driving at day & night, watching television etc. For activities required intermediate vision e.g. seeing car dashboard while driving, cooking & using computers, the satisfaction level was 100%, 100% and 96% respectively.



MODELS AVAILABLE

Model	Labeled Cylinder Power	Cylindrical Power		Recommended Range of Corneal Astigmatism Correction
		At IOL Plane	At Corneal Plane*	
PLHFD6T	CYL D: 1.00 D	1.00 D	0.68 D	0.25 D to 0.86.0 D
PLHFD6T	CYL D: 1.50 D	1.50 D	1.03 D	0.87 D to 1.25 D
PLHFD6T	CYL D: 2.25 D	2.25 D	1.54 D	1.26 D to 1.75 D
PLHFD6T	CYL D: 3.00 D	3.00 D	2.05 D	1.76 D to 2.25 D
PLHFD6T	CYL D: 3.75 D	3.75 D	2.57 D	2.26 D to 2.75 D
PLHFD6T	CYL D: 4.50 D	4.50 D	3.08 D	2.76 D to 3.25 D
PLHFD6T	CYL D: 5.25 D	5.25 D	3.60 D	3.26 D to 3.75 D
PLHFD6T	CYL D: 6.00 D	6.00 D	4.11 D	3.76 D and above

To choose suitable EYECRYL SERT TORIC model, please logon to



SPECIFICATION	P L H F D 6	P L H F D 6 T	
MATERIAL	Hydrophobic Acrylic Containing Natural Yellow Chromophore		
OPTIC TYPE	Single Piece 360° Square Edge with Aspheric Optic		
OPTIC SIZE	6.00 mm		
OVERALL SIZE	13.00 mm		
ANGULATION	0°		
ACD	5.28		
REFRACTIVE INDEX	1.524		
RECOMMENDED ULTRASOUND A-CONSTANT	SRK-T 119.00		
RECOMMENDED OPTICAL A-CONSTANTS	HAIGIS: α_0 :1.574, α_1 :0.400, α_2 :0.100	HOFFER - Q ACD: 5.78	
	Holl 1: 1.99	Holl 2: 2.02	SRK-T 119.40
	SRK-II 119.90	Barrett: 2.09	
DIOPTR RANGE	<p>+7.0 D to +30.0 D (with 0.5 D step) & 0.00 D to +6.5 D and +30.5 D to +40.0 D will be available on customization on case to case basis as per manufacturing possibility"</p>		
CYLINDER RANGE (For Toric Version)	1.0 D to 1.5 D (with 0.5 D step), 1.5 to 6.0 D (with 0.75 D step).		
	6.75 to 15.0 D (with 0.75 D step) customization as per manufacturing possibility		
IMPLANTATION SITE	Capsular Bag		
DELIVERY SYSTEM	Pre-loaded delivery system		
STERILIZATION	EO		

biotech

25 YEARS OF
innovation



Scan For Biotech Connect App

BIOTECH VISION CARE PVT. LTD.

Block 1, Abhisree Corporate Park, Opp. Swagat Bunglow BRTS Stop

Bopal - Ambli Road, Ahmedabad - 380 058 Gujarat India | Tel: +91 79 66823000

Email: Domestic: indiasales@biotechhealthcare.com | International: intlsales@biotechhealthcare.com | www.biotechhealthcare.com