

CE
2460

EYECRYL™ PLUS

Monofocal Hydrophilic Acrylic
Foldable Intraocular Lens



biotech

VISIONCARE

EYECRYL™ PLUS

Monofocal Hydrophilic Acrylic
Foldable Intraocular Lens


FEATURES

- 360 ° square edge
- 5 ° angulation ensures better contact with posterior capsule
- UV blocking Material
- Medical implantable grade hydrophilic material
- Packed in medical grade blister pack

SPECIFICATIONS

MATERIAL	Hydrophilic Acrylic CQ UV									
OPTIC TYPE	Single Piece, 360° Square Edge									
OPTIC SIZE	6.00 mm									
OVERALL SIZE	12.50 mm									
ANGULATION	5°									
ACD	4.96									
REFRACTIVE INDEX	1.462									
RECOMMENDED ULTRASOUND A-CONSTANT	SRK-T 118.0									
RECOMMENDED OPTICAL A-CONSTANTS	<table><tr><td>SRK - T 118.4</td><td>SRK - II 118.5</td><td>Holl 1 Const SF : 1.42</td></tr><tr><td colspan="2">HOFFER Q ACD : 5.19</td><td>HAIGIS a0:0.940, a1:0.40, a2:0.10</td></tr><tr><td>Barrett 1.57</td><td>Hill-RBF 118.3</td><td>Olsen ACD 4.20</td></tr></table>	SRK - T 118.4	SRK - II 118.5	Holl 1 Const SF : 1.42	HOFFER Q ACD : 5.19		HAIGIS a0:0.940, a1:0.40, a2:0.10	Barrett 1.57	Hill-RBF 118.3	Olsen ACD 4.20
SRK - T 118.4	SRK - II 118.5	Holl 1 Const SF : 1.42								
HOFFER Q ACD : 5.19		HAIGIS a0:0.940, a1:0.40, a2:0.10								
Barrett 1.57	Hill-RBF 118.3	Olsen ACD 4.20								
DIOPTR RANGE	-10.0 D to +30.0 D (with 0.5 D steps from +1.0D to +30.0D)									
IMPLANTATION SITE	Capsular Bag									
STERILIZATION	Steam									
SHELF LIFE	4 years from date of manufacture									

Model: 600



Scan For Biotech Connect App