


# BIO-HYALUR LVD

CE  
2460

Tranquil Surgeries. Calm & Clear Eyes.

Sodium Hyaluronate + Sulphate with Sulfate  
Viscoelastic Solution for Ophthalmic  
Anterior Segment Procedures



biotech

VISIONCARE

**Bio-Hyalur LVD** is used in Anterior Segment procedures including complicated Cataract Surgeries

## FEATURES

- + Protects the corneal endothelium and other ocular tissues
- + Maintains a deep chamber during anterior segment surgeries
- + Enhances visualization during the surgical procedure
- + The viscoelasticity of the solution maintains the normal position of the vitreous face, thus preventing formation of a postoperative flat chamber.

## TECHNICAL SPECIFICATIONS

<b>RAW MATERIAL</b>	Sodium Hyaluronate content (3.0% W/V) & Sodium Chondroitin Sulfate content (4.0% W/V)
<b>VISCOSITY (@ 25°C)</b>	40,000 ± 20,000 mPas
<b>pH</b>	6.5 to 7.6
<b>OSMOLARITY</b>	285 – 365 mOsm/kg
<b>SHELF LIFE</b>	2 years from the date of manufacturing
<b>STORAGE</b>	Store between 2°C to 8°C, Do not Freeze, Protect from light
<b>APPLICATION / INTENDED USE</b>	Anterior segment procedures

### BIOTECH VISION CARE PVT. LTD.

Block 1, Abhishree Corporate Park, Opp. Swagat Bungalow BRTS Stop  
Bopal - Ambli Road, Ahmedabad - 380 058 Gujarat India | Tel: +91 79 66823000

Email: Domestic: [indiasales@biotechhealthcare.com](mailto:indiasales@biotechhealthcare.com) | International: [intlsales@biotechhealthcare.com](mailto:intlsales@biotechhealthcare.com) | [www.biotechhealthcare.com](http://www.biotechhealthcare.com)