

GLODERM *Comfort*

Spirit of Natural Beauty

NEW GENERATION DERMAL FILLERS

GLODERM is a new and innovative range of smooth hyaluronic acid gels designed to enhance natural beauty.

Helps to volumize the face, remove sagging, reduce the appearance of wrinkles and shape the facial contours.




GLODERM *Comfort* THE SCIENCE BEHIND



Produced with **Golden Ratio Technology**



High molecular weight HA used in **GLODERM**, cross-linked with BDDE to obtain a three-dimensional gel with short and long chains in optimal ratio.



Monophasic, smooth colourless gel

12 L

12 L

A gel with high biocompatibility, easy to shape, homogeneous and desired viscosity is obtained.

HA concentration (mg/ml)	12
Level of cross-linking	●
Viscosity	Low
Suggested injection dept	Superficial
Pack	1 x 1 ml
Lidocain	0.3 %
Suggested indications	Smoker lines, crow's feet, neck lines, decolette lines, superficial frontal wrinkles
Needle size	30 G