

# GLODERM *Comfort*

## *Spirit of Natural Beauty*

### NEW GENERATION DERMAL FILLERS

GLODERM is a new and innovative range of smooth hyaluronic acid gels designed to enhance natural beauty.

Helps to volumize the face, remove sagging, reduce the appearance of wrinkles and shape the facial contours.



### GLODERM *Comfort*


## THE SCIENCE BEHIND



Produced with **Golden Ratio Technology**



High molecular weight HA used in **GLODERM**, cross-linked with BDDE to obtain a three-dimensional gel with short and long chains in optimal ratio.



Monophasic, smooth colourless gel

15L

15L

A gel with high biocompatibility, easy to shape, homogeneous and high viscosity is obtained.

HA concentration (mg/ml)	15
Level of cross-linking	●○○○
Viscosity	Low
Suggested injection depth	Superficial to mid dermis
Pack	1 x 1 ml
Lidocaine	0.3 %
Suggested indications	Defining philtrum and borderline of lips, wrinkles of forehead, under brows, hand rejuvenation
Needle size	30G