

# GLODERM *Comfort*

## *Spirit of Natural Beauty*

### NEW GENERATION **DERMAL FILLERS**

GLODERM is a new and innovative range of smooth hyaluronic acid gels designed to enhance natural beauty.

Helps to volumize the face, remove sagging, reduce the appearance of wrinkles and shape the facial contours.



# GLODERM *Comfort*


## THE SCIENCE BEHIND



Produced with **Golden Ratio Technology**



High molecular weight HA used in **GLODERM**, cross-linked with BDDE to obtain a three-dimensional gel with short and long chains in optimal ratio.



24L

24L

A gel with high biocompatibility, easy to shape, homogeneous and high viscosity is obtained.

|                          |   |
|--------------------------|---|
| HA concentration (mg/ml) | 20  |
| Level of cross-linking   | ● ●   |
| Viscosity                | Mid   |
| Suggested injection dept | Subcutaneous  |
| Pack                     | 1 x 1 ml  |
| Lidocain                 | 0.3 %   |
| Suggested indications    | Lip augmentation, tear trough<br>moderate nasolabial folds,<br>marionette lines |
| Needle size              | 27G   |