

GLODERM

Spirit of Natural Beauty

NEW GENERATION **DERMAL FILLERS**

GLODERM is a new and innovative range of smooth hyaluronic acid gels designed to enhance natural beauty.

Helps to volumize the face, remove sagging, reduce the appearance of wrinkles and shape the facial contours.



GLODERM FILLERS


THE SCIENCE BEHIND



Produced with **Golden Ratio Technology**



High molecular weight HA used in **GLODERM**, cross-linked with BDDE to obtain a three-dimensional gel with short and long chains in optimal ratio.



20L

20L

A gel with high biocompatibility, easy to shape, homogeneous and high viscosity is obtained.

HA concentration (mg/ml)	20
Level of cross-linking	●●○
Viscosity	mid
Suggested injection depth	Mid dermis
Lifting capacity	low
Pack	1 x 1 ml
Suggested indications	Fine lines, nasolabial folds, marionette lines, lip contour, under eye pits
Needle size	27G ½"